



# Ligature Resistant Products

supplying products to create a safer environment

Call: 03300 536020

Email: [info@anti-ligature-shop.co.uk](mailto:info@anti-ligature-shop.co.uk)

Visit: [www.anti-ligature-shop.co.uk](http://www.anti-ligature-shop.co.uk)

## Fixing Reduced Ligature Cutter Box (2 Knife Version)

By Anti-Ligature Shop

**Product Code: ALSLCEB**

Using existing holes, fix to the wall with suitable wall fixings.

Knives are held in position using magnets, position of magnets is shown.

W 300 x L 370 x D 50mm



Place the knives on top of the magnets to hold them in position.

### Please note:

- Ensure cabinet is locked at all times.
- Knives should be replaced after use.
- On completion of installation it is recommended that a small bead of Anti-Pick 110 Sealant is used around the outer frame.
- **Do not** use solvent based products on this product.



### Cleaning:

Wipe down using anti-bacterial products.

### NBS Specification

The anti-ligature design features of this product are not intended to replace or substitute the need for necessary supervision of those who may be at risk, or for other necessary protective measures to be taken in the specific circumstances of usage. Anti-Ligature Shop does not claim any product to be totally ligature free or make claim that they cannot be damaged with regard to anti-vandal feature. It is therefore the customers responsibility to ensure the suitability of the product for the environment in which they are to be fitted.

Call: 03300 536020

Email: [info@anti-ligature-shop.co.uk](mailto:info@anti-ligature-shop.co.uk)

Visit: [www.anti-ligature-shop.co.uk](http://www.anti-ligature-shop.co.uk)



51 Haydock street | Soteria House | Bewsey & Whitecross  
| Warrington | Cheshire | WA2 7UW